



Zulassungen / Konformitäten

EPDM (peroxidisch vernetzt)

FDA konform (§ 21 CFR 177.2600)
USP Class VI
3 A Sanitary Standard, Number 18-03 (Class II)
DVGW (KTW, W270 und W534)
KTW
EN 681-1 WB

FKM/FEP ummantelt

FDA konform (§ 21 CFR 177.1550)
USP Class VI
3 A Sanitary Standard, Number 18-03 (Class I)

weitere Werkstoffe auf Anfrage

Verpackung

einzeln verpackt

Approvals / Conformities

EPDM (crosslinked peroxidic)

FDA compliant (§ 21 CFR 177.2600)
USP Class VI
3 A Sanitary Standard, Number 18-03 (Class II)
DVGW (KTW, W270 and W534)
KTW
EN 681-1 WB

FKM/FEP encapsulated

FDA compliant (§ 21 CFR 177.1550)
USP Class VI
3 A Sanitary Standard, Number 18-03 (Class I)

further materials on request

Packaging

wrapped separately



weitere Abmessungen auf Anfrage | further dimensions on request

DIN	Abmessung dimension						EPDM	FKM/ FEP **		
DN010	13.0 x 1.5						703297	705087		
DN015	19.0 x 1.5						703298	704162		
DN020	23.0 x 1.5						703299	705088		
DN025	29.0 x 1.5						703300	704160		
DN032	35.0 x 1.5						703301	705089		
DN040	41.0 x 1.5						703302	704161		
DN050	53.0 x 1.5						703303	704521		
DN065	70.0 x 2.0						703304	704638		
DN080	85.0 x 2.0						703305	705090		
DN100	104.0 x 2.0						703306	704504		
DN125	129.0 x 2.0						705664	704605		
DN150	154.0 x 2.0						704926	704606		
DN200*	204.0 x 2.0						708099	705495		
DN250*	254.0 x 2.0						706684	712119		
DN300*	306.0 x 3.0						708100	712042		

ISO	Abmessung dimension						EPDM	FKM/ FEP **		
DN008	13.5 x 1.6						704739	705091		
DN010	17.2 x 1.6						704738	705092		
DN015	21.3 x 1.6						703307	705093		
DN020	26.9 x 1.6						703308	705094		
DN025	33.7 x 2.0						703309	705095		
DN032	42.4 x 2.0						703310	705096		
DN040	48.3 x 2.0						703311	704884		
DN050	60.3 x 2.0						703312	704522		
DN065	76.1 x 2.0						703313	705097		
DN080	88.9 x 2.3						703314	705098		
DN100	114.3 x 2.3						703887	701497		
DN125*	139.7 x 2.6						703241	706737		
DN150*	168.3 x 2.6						703242	703789		
DN200*	219.1 x 2.6						703791	706738		
DN250*	273.0 x 3.0						703790	706739		
DN300*	323.9 x 3.0						703792	712043		

ASME BPE	Abmessung dimension						EPDM	FKM/ FEP **		
DN 1/2"	12.7 x 1.65						704737	705099		
DN 3/4"	19.05 x 1.65						704736	705100		
DN 1"	25.4 x 1.65						703315	705101		
DN 1 1/2"	38.1 x 1.65						703316	705102		
DN 2"	50.8 x 1.65						703317	705103		
DN 2 1/2"	63.5 x 1.65						703318	705104		
DN 3"	76.2 x 1.65						703319	705105		
DN 4"	101.6 x 2.11						703320	705106		

*) nicht in der DIN-Norm enthalten | non-standard

**) ummantelt | encapsulated